

### In vitro ADME & PK

# $pK_a$ and log P

## Background Information



<sup>6</sup>PK<sub>a</sub> affects solubility, permeability, log D and oral absorption by modulating the distribution of neutral and charged species.<sup>2</sup>

<sup>1</sup>Di L and Kerns EH. (2003) *Current Opinion in Chemical Biology* **7**; 402-408.

- The pK<sub>a</sub> is the pH at which the molecule is 50% protonated.
- Log P (or partition co-efficient) is a measure of the lipophilicity of a compound.
- Cyprotex's pK<sub>a</sub> and log P determination uses UV-metric and pH-metric technology developed by Sirius, which is considered to be a 'gold standard' method for determining these properties.
- In UV-metric methods, the pK<sub>a</sub> is measured by analysing changes in multi-wavelength UV spectra during acid-based titration of the sample. UV-metric pK<sub>a</sub> methods work for compounds with pH-sensitive chromophores.
- In pH-metric methods, pK<sub>a</sub> is measured by titrating a solution of the sample in water or solvent with acid and base, and calculating the pK<sub>a</sub> from the shape of the titration. pH-metric methods work for any ionisable compound, but require more sample than UV-metric methods.
- The pH-metric method is also used to measure log P in a two-phase acid-base titration in the presence of octanol.

#### Protocol

#### Method

Fast UV titration for  $pK_a$  UV-metric titration for  $pK_a$  Potentiometric (pH-metric) titration for  $pK_a$  and log P

#### Instrument SiriusT3

#### Test Article Requirements

3-5 μL of 10 mM stock solution (UV-metric) 1 mg solid compound (pH-metric)

Partition Solvent used for Log P Determination n-Octanol (others available on request)

#### **Data Delivery**

pK<sub>a</sub> log P (optional) Standard error RMSD Calculated log D at pH<sub>7.4</sub> (based on pK<sub>a</sub> and log P)

#### log P is determined from the shape of titration curves obtained in dual-phase titrations.



#### pK<sub>and</sub> log P

pK<sub>a</sub> measurements are determined using the SiriusT3 instrument from Sirius-Analytical using either a UV-metric or pH-metric approach.

#### Figure 1

UV-metric method for measuring pK<sub>a</sub> values.



UV-metric methods provide  $pK_a$  results for samples with chromophores whose UV absorbance changes as a function of pH.

On SiriusT3, the Fast UV method measures absorbance at 250 wavelengths and 54 pH values in a buffered solution in about 5 minutes. The slower UV-metric method in unbuffered solution extends the pH range below 1 or above 13. The 3D graph shows data from the measurement of labetalol pK<sub>a</sub>s. The other graphs are 2D projections showing change in absorbance vs. pH and vs. wavelength, with percent species and molar absorbance coefficients overlain.





#### Figure 2

pH-metric method for measuring pK<sub>a</sub> values.



pH-metric methods are based on potentiometric acid-base titration. Results are obtained by a complex computational process. The pH of each point in the titration curve is calculated using equations that contain  $pK_{a'}$  and the calculated points are fitted to the measured curve by manipulating the  $pK_a$  value. The  $pK_a$  that provides the best fit is taken to be the measured  $pK_a$ . pH-metric methods will measure all  $pK_a$  between 2 and 12, provided the sample is in solution throughout the experiment.

#### Figure 3

pH-metric method for measuring log P values



In the pH-metric method for log P, a weighed sample is dissolved in a two-phase water-octanol system, and titrated over a pH range (typically 2 to 12 for bases and ampholytes, 12 to 2 for acids). Although the solution becomes opaque during stirring, the pH electrode continues to measure pH of the aqueous component of the solution. Results are obtained by a complex computational process. The pH of each point in the titration curve is calculated using equations that contain pK<sub>a</sub> and P, and the calculated points are fitted to the measured curve by manipulating the P value. The P that provides the best fit is taken to be the measured P value, which is reported as the logarithm, i.e. log P. As well as log P, the log D value as a function of pH is determined from the data.

#### References

<sup>1</sup> Di L and Kerns EH. (2003) Current Opinion in Chemical Biology 7; 402-408.