

Validation of a serum-free approach to facilitate the development of high-throughput immunogenicity screening assays *in vitro*.

Stephen Madden, Emma Shardlow, James Copley & Paul Walker

Introduction

- The generation of anti-drug antibodies (ADA) towards a therapeutic agent can have severe implications for both drug safety and efficacy
- Classical *in vivo* models often fail to fully recapitulate the complexity of the human immune response and thus may not adequately predict the immunogenic potential of candidates entering the pipeline¹
- As part of an early de-risking strategy, we present a high-throughput, *in vitro* approach to immunogenicity screening using two distinct models: i) peripheral blood mononuclear cells (PBMC) and ii) monocyte-derived dendritic cells (moDCs) isolated from the blood of healthy human donors
- In order to standardise our platform, we adopted a serum-free approach to cell culture which can be used to assess the immunogenicity of myriad modalities including biologics, oligonucleotides and small molecules

Method

- To investigate the impact of serum removal upon cellular viability, proliferation and cytokine release, PBMC and moDCs were cultured using a commercially available serum-free media (X-Vivo 15, Lonza) or RPMI 1640 supplemented with 10% human AB serum (serum-containing (SC) conditions)
- CD14⁺ monocytes were seeded in 384-well plates and differentiated into moDCs using IL-4 & GM-CSF for 6 days prior to compound exposure. On day 6, moDCs were dosed with various concentrations of test articles in order to produce an 8-point dose response curve
- Following a 24hr incubation period, supernatants were collected via centrifugation and cells were assayed for viability (ATP content) using a CellTiter-Glo® luminescent assay (Promega). Cytokine analysis was performed upon collected supernatants using a multiplex cytokine bead array (BD Bioscience) according to the manufacturers instructions and data collected using an Intellicyt iQue 3 flow cytometer (Sartorius)
- For PBMC, cells were seeded in 384-well plates and exposed to compounds after a 24hr resting period. Lymphocytic proliferation and cell viability were measured in tandem after 72hrs using a resazurin based fluorescent assay (AlamarBlue™, Thermo Fisher scientific) and CellTiter-Glo® kit. Supernatants were collected and analysed for cytokine release after 72hrs as previously described

Serum-free (SF) conditions improve the viability and morphology of moDCs in culture




Figure 1: Brightfield images (40X mag.) illustrating the morphological differences between cells cultured in SC conditions (Figs. A-C) versus those cultured using SF conditions (Figs. D-F). Figs. A & D represent moDCs exposed to vehicle only (HBSS) while B & E and C & F highlight the typical phenotypic shift observed upon treatment with LPS (1µg/mL) and ODN 2395 (10µM) respectively for 24hrs. Evidence of cellular debris was reduced under SF conditions (Figs. A vs. D) and this was reflected by a moderate increase in ATP content (data not shown).

SF conditions enhance levels of cytokine release associated with PBMC proliferation while increasing the sensitivity of the plate-based lymphocytic proliferation assay




Figure 2: PBMC were exposed to PHA (0.01-25µg/mL) for 72hrs under SC (Figs. A-C) and SF conditions (Figs. D-F). Lymphocytic proliferation was assessed by measuring resazurin reduction (Figs B & E) and cellular ATP content (Figs A & D). Cytokine secretion was measured by flow cytometry using a multiplex CBA kit and the levels of IL-2, IFN-γ, IL-6 and TNF after exposure to 2.5 µg/mL PHA are shown (Figs. C & F). Red dotted lines represent both the lower solvent limit and an upper limit that conforms to the standard currently applied in assays used to evaluate compound immunogenicity (stimulation index (SI) - 2). Error bars show ±SD of triplicate intra assay replicates (n=1 donor). Unfilled circles indicate data points that were excluded from the final analysis.

- MEC values for lymphocytic proliferation generated under SF conditions were lower than those obtained from SC cultures (AlamarBlue™ – 0.21 & 0.33 µg/mL respectively, ATP content – 0.19 & 0.29 µg/mL respectively)
- MEC values for IL-2 release were reduced in the absence of serum (0.25µg/mL vs 0.42 µg/mL) and levels were ca 60 fold greater under SF conditions

Measuring PBMC cytokine release and proliferation under SF conditions can be used to evaluate the immunogenic potential of various drug modalities *in vitro*

	Immunogenic <i>in vivo</i> ?	IL-6	IL-2	IFN-γ	Proliferation
Biologics					
1. PHA	✓	↑ <0.004µg/mL	↑ <0.004µg/mL	↑ <0.004µg/mL	↑ <0.004µg/mL
2. Bevacizumab	ADA rate – 0% ³	NR	NR	NR	NR
3. Adalimumab	ADA rate – 28% ³	NR	NR	NR	NR
Oligonucleotides					
4. ODN2395 (Class C)	✓ - IL-6 secretion ¹	↑ 0.009µM	NR	NR	↑ 0.5µM
Small molecules					
5. CL075	✓	↑ 0.29µg/mL	↑ 0.62µg/mL	NR	↑ 3.1µg/mL

Table 1: MEC values derived from lymphocytic proliferation (AlamarBlue™) and cytokine release data for various modalities including biologics (1-3), oligonucleotides (4) and small molecules (5). IL-2 and IFN-γ were selected as key markers of T cell activation while IL-6 was used as an indicator of a pro-inflammatory response². PHA is a T cell mitogen and is known to stimulate the production of large quantities of IL-2 and IFN-γ. Bevacizumab and adalimumab are antibodies that target VEGF and TNF-α respectively and are generally considered to present a low-medium risk in terms of the generation of ADA in a clinical setting³. ODN2395 is an oligonucleotide that is known to stimulate significant IL-6 production via interaction with B cell and DC targets *in vivo*¹. CL075 is a TLR8 agonist that activates the NF-κB pathway culminating in the production of pro-inflammatory signals such as IL-6.⁴ NR refers to a compound that is non-responsive in that feature of the assay and the arrow indicates the direction of the response.

SF moDCs can be used to evaluate both the toxicity and propensity of a compound to elicit a pro-inflammatory response at an early stage

	Immunogenic <i>in vivo</i> ?	IL-6	TNF	IL-10	Cellular ATP
Biologics					
1. LPS EK	✓	↑ <0.0004µg/mL	↑ 0.001µg/mL	↑ 0.004µg/mL	NR
2. MDP	✓	↑ 0.01µg/mL	↑ 0.14µg/mL	NR	NR
Oligonucleotides					
3. ODN2395 (Class C)	✓	↑ 0.85µM	↑ 0.91µM	↑ 5.95µM	NR
4. ODN1668 (Class B)	✓	↑ 2.29µM	↑ 1.40µM	NR	NR
5. ssRNA40/LyoVec™	✓	↑ 0.13µg/mL	↑ 0.33µg/mL	NR	NR
Small molecules					
6. Resiquimod (R848)	✓	↑ 0.05µg/mL	↑ 0.29µg/mL	↑ 1.43µg/mL	NR
7. CL075	✓	↑ 0.06µg/mL	↑ 0.05µg/mL	↑ 0.32µg/mL	NR
8. Pindolol	*	NR	NR	NR	NR
9. Lapatinib	*	NR	NR	NR	↓ 4.65µM

Table 2: MEC values derived from cellular ATP and cytokine release data for various modalities including biologics (1-2), oligonucleotides (3-5) and small molecules (6-9). Compounds 1-7 are all known to be immunogenic whereas compound 8 has no known propensity to elicit an immune response *in vivo*. While lapatinib exhibits toxicity towards moDCs due to its chemotherapeutic nature, it has also been associated with delayed-onset hypersensitivity reactions; however, these are idiosyncratic and immune activation only occurs in susceptible individuals. NR refers to a compound that is non-responsive in that feature of the assay and the arrow indicates the direction of the response.

Conclusions

- Serum-free conditions provide a suitable environment in which to culture primary immune cells intended for downstream application in immunogenicity assays
- Eliminating serum enhances cytokine signal and assay performance without negatively impacting cell viability
- Measuring viability, lymphocytic proliferation and cytokine release in tandem provides the greatest mechanistic insight into the immunogenic affects elicited by candidate drugs
- Both moDCs and PBMCs can be used to accurately predict the likelihood of a pro-inflammatory response in a rapid, high-throughput format
- PBMC show greater sensitivity with regards to the prediction of a pro-inflammatory response where multiple cell types are involved in directing an immune response i.e. ODN2395

References

1. Burel *et al.*, (2022) Specific Inflammatory Responses in Human Volunteer Peripheral Blood Mononuclear Cells. *Nucleic Acid Ther.*, **32(6)**, 457-472.
2. Joubert *et al.*, (2016) Use of *In Vitro* Assays to Assess Immunogenicity Risk of Antibody-Based Biotherapeutics. *PLoS One*, **11(8)**, e0159328.
3. Vaisman-Mentesh *et al.*, (2020) The Molecular Mechanisms That Underlie the Immune Biology of Anti-drug Antibody Formation Following Treatment With Monoclonal Antibodies. *Front. Immunol.*, **11**, 1951.
4. Spranger *et al.*, (2010) Generation of Th1-polarizing dendritic cells using the TLR7/8 agonist CL075. *J. Immunol.* **185**, 738-747.